

The Engineering of the Walk'n'Chair

Mike
Johnson M.D.


Why does the
Walk'n'Chair morph into
so many configurations?

- Because people have different mobility limitations.
- Because people have different environments to contend with.

**Outdoor and indoor
environments
require an adaptive
mobility device.**

And some serious
engineering.

The Walk'n'Chair Journal



The (Adaptive) Outdoor Activity
Extender







Walk'n'Chair™

by
BROOKFIELD HUNTER INCORPORATED

The outdoor environment is not always conducive to mobility aids.

As the outdoor terrain becomes more irregular, it is second nature to take more care and go slowly.

Theoretically, increasing the rolling resistance of the wheels and increasing the turning radius as terrain gets more irregular mimics what people do when the going gets rough. They slow down and take wider turns.

The configurations of Walk'n' Chair.				
Purpose	Rough Terrain Walker	Smooth Terrain Walker	Rough Terrain Wheel Chair	Smooth Terrain Wheel Chair with foot rest
Turning	Wide radius	0 radius	Wide radius	Narrow radius
Rolling Resistance	+++ Resistance	+ Resistance	+++ Resistance	+ Resistance

Rolling Resistance



Rolling resistance is related to the tire material, wheel radius and load on the wheel base (Lin et al., 2015; Mediola et al., 2014; Sprigle and Huang 2015). Most of the resistance to rolling is due to “tire scrub” and “caster shimmy” not bearing resistance².

How Walk'n'Chair changes Rolling Resistance

Changing the configuration of the Walk'n'Chair alters the weight distribution of the load relative to the front and back wheels. The large wheels have much less rolling friction than to the smaller caster wheels (due to size, materials, caster shimmy).

When we change the Walk'n'Chair from a **Smooth Terrain Walker** to a **Rough Terrain Walker**, the center of gravity moves from the large wheels (low resistance) to the casters (high resistance)

This increases the rolling resistance to match the caution one normally uses in rough terrain. Increased resistance for rough terrain gives the Walk'n'Chair a solid feel.

	
<p>Rough Terrain Walker</p>	<p>Smooth Terrain Walker</p>
<p>Center of gravity is closest to high rolling resistance caster</p>	<p>Center of gravity is closest to low resistance larger wheel</p>
<p>Resistance +++</p>	<p>Resistance +</p>

The Rough Terrain Walker increases the rolling resistance to match the slow and cautious strategy one assumes with rough ground.

Turning Radius (Maneuverability)

The Walk'n'Chair turning radius is increased for the **Rough Terrain Walker** configuration. This in

combination with increased rolling resistance makes the **Rough Terrain Walker** configuration a sturdy match for irregular terrain.

Effects of Wheel Camber

The inclination of the large rear wheels (camber) is usually reserved for athletics and wheelchair sports (Tsai et al., 2012). Camber increases lateral stability; the effects are significant from 0-15°. Camber allows increased turning velocities. Camber from 8 to 15° increases

rolling resistance (Tsai et al., 2012). A camber of 6° was found to be optimal in terms of lateral stability on smooth terrain (Mediola et al., 2014).

The Walk'n'Chair large wheel camber is 12°. This enhances stability especially when moving over lateral slopes.

The Walk'n'Chair was designed to traverse the obstacles in your way.



Rough Terrain Walker

Large wheel traverses a curb obstacle. Curb can be $\frac{1}{2}$ of radius of wheel.



Smooth Terrain Walker

Small caster wheel traverses a curb obstacle.

Obstacle Traverse

The largest obstacle height safely traversed is related to the wheel radius. A curb half of the radius of the forward wheel is easily and safely crossed.

The ***Rough Terrain Walker*** and ***Rough Terrain Wheelchair*** configurations allow the traverse of rough terrain and curbs.

Conclusion

Several principles of wheel chair design were extracted from recent engineering literature¹⁻⁸. These principles were applied to the Walk'n'Chair.

The Walk'n'Chair modifies the center of gravity to match rolling resistance to the terrain. Turning radius is modified to optimize for smooth indoor or rough outdoor terrain. Stability is increased by reducing turning radius and adding camber to the rear wheels. Rough terrain features are

traversed with larger wheels.

References

1. Caspall JJ, Selgsohn E, Dao PV, Springle S (2013) Changes in inertia and effect on turning effort across different wheelchair configurations. J Rehabil Res Dev 50(10):1353-62.
2. Lin JT, Huang M, Springle S (2015) Evaluation of wheelchair resistive forces during straight and turning trajectories across different wheelchair configurations using free-wheeling coast-down test. J Rehab Research Devel 52 (7): 763-774.
3. Mediola FO, Elui VMC, Santana CDS, Fortulan CA (2014) Aspects of manual wheelchair configuration affecting mobility: A review. J Phys Ther Sci 26:313-318.
4. Sauret C, Bascou J, de Saint Rémy N, Pillet H, Vaslin P, Lavaste F. (2012) Assessment of field rolling resistance of manual wheelchairs. J Rehabil Res Dev. 2012;49(1):63-74.
5. Springle S, Huang M (2015) Impact of mass and weight distribution on manual wheelchair propulsion torque. Assist Technol 27(4) 226-235.
6. Thomas L, Borisoff J, Sparrey CJ (2017) Quantifying the effects of on-the-fly changes of seating configuration on the stability of a manual wheelchair. Conf Proc IEEE Eng Biol Soc 2017: 1897-1900.
7. Tsai CY, Lin CJ, Huang YC, Lin PC, Su FC (2012) The effects of rear-wheel

camber on the kinematics of upper extremity during wheelchair propulsion. Biomed Eng Online. 2012 Nov 22;11:87.

8. Zepeda R, Chan F, Sawatzky B (2017) The effect of caster wheel diameter and mass distribution on drag forces in manual wheelchairs. J Rehabil Res Dev 53(6) 893-900

Check out the Walk'n'Chair at www.walknchair.com

Or email us: USA sales@walknchair.com
 Canada sales@walknchair.com

Or call us: USA 1(406) 925 5663
 Canada 250 324 1285